

Product datasheet for **TL300558**

VPS26B Human shRNA Plasmid Kit (Locus ID 112936)

Product data:

Product Type:	shRNA Plasmids
Product Name:	VPS26B Human shRNA Plasmid Kit (Locus ID 112936)
Locus ID:	112936
Synonyms:	Pep8b
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	VPS26B - Human, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 112936). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	NM_052875 , NM_052875.1 , NM_052875.2 , NM_052875.3 , NM_052875.4 , BC098386 , BC098386.1 , BC007799 , BC009747 , BC013255 , BC014128 , NM_052875.5
UniProt ID:	Q4G0F5
Summary:	Acts as component of the retromer cargo-selective complex (CSC). The CSC is believed to be the core functional component of retromer or respective retromer complex variants acting to prevent missorting of selected transmembrane cargo proteins into the lysosomal degradation pathway. The recruitment of the CSC to the endosomal membrane involves RAB7A and SNX3. The SNX-BAR retromer mediates retrograde transport of cargo proteins from endosomes to the trans-Golgi network (TGN) and is involved in endosome-to-plasma membrane transport for cargo protein recycling. The SNX3-retromer mediates the retrograde transport of WLS distinct from the SNX-BAR retromer pathway. The SNX27-retromer is believed to be involved in endosome-to-plasma membrane trafficking and recycling of a broad spectrum of cargo proteins. The CSC seems to act as recruitment hub for other proteins, such as the WASH complex and TBC1D5. May be involved in retrograde transport of SORT1 but not of IGF2R. Acts redundantly with VSP26A in SNX-27 mediated endocytic recycling of SLC2A1/GLUT1 (By similarity).[UniProtKB/Swiss-Prot Function]



[View online »](#)


shRNA Design:

These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com. If you need a special design or shRNA sequence, please utilize our [custom shRNA service](#).


Performance Guaranteed:

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.


For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).

Product images:


GFP signal was observed under microscope at 48 hours after transduction of TL300558A virus into HEK293 cells. TL300558A virus was prepared using lenti-shRNA TL300558A and [TR30037] packaging kit.



GFP signal was observed under microscope at 48 hours after transduction of TL300558B virus into HEK293 cells. TL300558B virus was prepared using lenti-shRNA TL300558B and [TR30037] packaging kit.



GFP signal was observed under microscope at 48 hours after transduction of [TL300558C] virus into HEK293 cells. [TL300558C] virus was prepared using lenti-shRNA [TL300558C] and [TR30037] packaging kit.



GFP signal was observed under microscope at 48 hours after transduction of [TL300558D] virus into HEK293 cells. [TL300558D] virus was prepared using lenti-shRNA [TL300558D] and [TR30037] packaging kit.