

## Product datasheet for **TL300570**

### Bestrophin 3 (BEST3) Human shRNA Plasmid Kit (Locus ID 144453)

#### Product data:

Product Type:	shRNA Plasmids
Product Name:	Bestrophin 3 (BEST3) Human shRNA Plasmid Kit (Locus ID 144453)
Locus ID:	144453
Synonyms:	VMD2L3
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	BEST3 - Human, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 144453). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	<a href="#">NM_001282613</a> , <a href="#">NM_001282614</a> , <a href="#">NM_001282615</a> , <a href="#">NM_001282616</a> , <a href="#">NM_032735</a> , <a href="#">NM_152439</a> , <a href="#">NR_104220</a> , <a href="#">NM_032735.1</a> , <a href="#">NM_032735.2</a> , <a href="#">NM_152439.1</a> , <a href="#">NM_152439.2</a> , <a href="#">NM_152439.3</a> , <a href="#">NM_001282615.1</a> , <a href="#">NM_001282616.1</a> , <a href="#">NM_001282614.1</a> , <a href="#">NM_001282613.1</a> , <a href="#">BC006440</a> , <a href="#">BC028087</a> , <a href="#">BC151138</a> , <a href="#">BC166663</a> , <a href="#">NM_001282615.2</a> , <a href="#">NM_001282616.2</a> , <a href="#">NM_032735.3</a>
UniProt ID:	<a href="#">Q8N1M1</a>
Summary:	BEST3 belongs to the bestrophin family of anion channels, which includes BEST1 (MIM 607854), the gene mutant in vitelliform macular dystrophy (VMD; MIM 153700), and 2 other BEST1-like genes, BEST2 (MIM 607335) and BEST4 (MIM 607336). Bestrophins are transmembrane (TM) proteins that share a homology region containing a high content of aromatic residues, including an invariant arg-phe-pro (RFP) motif. The bestrophin genes share a conserved gene structure, with almost identical sizes of the 8 RFP-TM domain-encoding exons and highly conserved exon-intron boundaries. Each of the 4 bestrophin genes has a unique 3-prime end of variable length (Stohr et al., 2002 [PubMed 12032738]; Tsunenari et al., 2003 [PubMed 12907679]).[supplied by OMIM, Mar 2008]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <a href="mailto:techsupport@origene.com">techsupport@origene.com</a> . If you need a special design or shRNA sequence, please utilize our <a href="#">custom shRNA service</a> .



[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).