

Product datasheet for **TL300582**

VEZT Human shRNA Plasmid Kit (Locus ID 55591)

Product data:

Product Type:	shRNA Plasmids
Product Name:	VEZT Human shRNA Plasmid Kit (Locus ID 55591)
Locus ID:	55591
Synonyms:	VEZATIN
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	VEZT - Human, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 55591). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	NM_017599 , NR_038241 , NR_038242 , NM_001352088 , NM_001352089 , NM_001352090 , NM_001352091 , NM_001352092 , NM_001352093 , NM_001352094 , NM_001352095 , NM_001352096 , NM_001352097 , NM_001352098 , NM_001352099 , NM_001352100 , NM_001352101 , NM_001352102 , NM_001352103 , NM_001352104 , NM_001352105 , NM_001352106 , NM_001352107 , NM_001352108 , NM_001352109 , NM_001352110 , NM_001352111 , NM_001352112 , NM_001352113 , NM_001352114 , NM_001352115 , NM_001352116 , NM_001352117 , NM_001352118 , NM_001352119 , NR_147917 , NR_147918 , NR_147919 , NR_147920 , NR_147921 , NR_147922 , NR_147923 , NR_147924 , NR_147925 , NR_147926 , NM_017599.1 , NM_017599.2 , NM_017599.3 , BC014443 , BC028597 , BC040845 , BC064939 , NM_017599.4
UniProt ID:	Q9HBM0
Summary:	This gene encodes a transmembrane protein which has been localized to adherens junctions and shown to bind to myosin VIIA. Examination of expression of this gene in gastric cancer tissues have shown that expression is decreased which appears to be related to hypermethylation of the promoter. Expression of this gene may also be inhibited by binding of a specific microRNA to a target sequence in the 3' UTR of the transcripts. A pseudogene of this gene is located on the X chromosome. Alternative splicing results in multiple transcript variants. [provided by RefSeq, May 2011]



[View online »](#)

- shRNA Design:** These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com. If you need a special design or shRNA sequence, please utilize our [custom shRNA service](#).
- Performance Guaranteed:** OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.
- For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).