

## Product datasheet for **TL300788**

### ZNF537 (TSHZ3) Human shRNA Plasmid Kit (Locus ID 57616)

#### Product data:

Product Type:	shRNA Plasmids
Product Name:	ZNF537 (TSHZ3) Human shRNA Plasmid Kit (Locus ID 57616)
Locus ID:	57616
Synonyms:	TSH3; ZNF537
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	TSHZ3 - Human, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 57616). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	<a href="#">NM_020856</a> , <a href="#">NR_138034</a> , <a href="#">NR_138035</a> , <a href="#">NR_138036</a> , <a href="#">NM_020856.1</a> , <a href="#">NM_020856.2</a> , <a href="#">NM_020856.3</a> , <a href="#">BC127095</a> , <a href="#">BC127096</a> , <a href="#">NM_020856.4</a>
UniProt ID:	<a href="#">Q63HK5</a>
Summary:	This gene encodes a zinc-finger transcription factor that regulates smooth muscle cell differentiation in the developing urinary tract. Consistent with this role, mice in which this gene has been inactivated exhibit abnormal gene expression in urinary tract smooth muscle cell precursors and kidney defects including hydronephrosis. The encoded transcription factor comprises a gene silencing complex that inhibits caspase expression. Reduced expression of this gene and consequent caspase upregulation may be correlated with progression of Alzheimer's disease in human patients. [provided by RefSeq, Jul 2016]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <a href="mailto:techsupport@origene.com">techsupport@origene.com</a> . If you need a special design or shRNA sequence, please utilize our <a href="#">custom shRNA service</a> .



[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).