

Product datasheet for **TL300807V**

VR1 (TRPV1) Human shRNA Lentiviral Particle (Locus ID 7442)

Product data:

Product Type:	shRNA Lentiviral Particles
Product Name:	VR1 (TRPV1) Human shRNA Lentiviral Particle (Locus ID 7442)
Locus ID:	7442
Synonyms:	VR1
Vector:	pGFP-C-shLenti (TR30023)
Format:	Lentiviral particles
Components:	TRPV1 - Human shRNA lentiviral particles (4 unique 29mer target-specific shRNA, 1 scramble control), 0.5 ml each, >10 ⁷ TU/ml.
RefSeq:	NM_018727 , NM_080704 , NM_080705 , NM_080706 , NM_080704.1 , NM_080704.2 , NM_080704.3 , NM_080705.1 , NM_080705.2 , NM_080705.3 , NM_080706.1 , NM_080706.2 , NM_080706.3 , NM_018727.1 , NM_018727.2 , NM_018727.3 , NM_018727.4 , BC132820 , BC136633 , NM_080704.4 , NM_080705.4
UniProt ID:	Q8NER1
Summary:	Capsaicin, the main pungent ingredient in hot chili peppers, elicits a sensation of burning pain by selectively activating sensory neurons that convey information about noxious stimuli to the central nervous system. The protein encoded by this gene is a receptor for capsaicin and is a non-selective cation channel that is structurally related to members of the TRP family of ion channels. This receptor is also activated by increases in temperature in the noxious range, suggesting that it functions as a transducer of painful thermal stimuli in vivo. Four transcript variants encoding the same protein, but with different 5' UTR sequence, have been described for this gene. [provided by RefSeq, Jul 2008]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com . If you need a special design or shRNA sequence, please utilize our custom shRNA service .



[View online »](#)

**Performance
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).