

## Product datasheet for **TL300883V**

### Two pore calcium channel protein 2 (TPCN2) Human shRNA Lentiviral Particle (Locus ID 219931)

#### Product data:

Product Type:	shRNA Lentiviral Particles
Product Name:	Two pore calcium channel protein 2 (TPCN2) Human shRNA Lentiviral Particle (Locus ID 219931)
Locus ID:	219931
Synonyms:	SHEP10; TPC2
Vector:	pGFP-C-shLenti (TR30023)
Format:	Lentiviral particles
Components:	TPCN2 - Human shRNA lentiviral particles (4 unique 29mer target-specific shRNA, 1 scramble control), 0.5 ml each, >10 <sup>7</sup> TU/ml.
RefSeq:	<a href="#">NM_139075</a> , <a href="#">NM_139075.1</a> , <a href="#">NM_139075.2</a> , <a href="#">NM_139075.3</a> , <a href="#">BC063008</a> , <a href="#">BC063008.1</a>
UniProt ID:	<a href="#">Q8NHX9</a>
Summary:	This gene encodes a putative cation-selective ion channel with two repeats of a six-transmembrane-domain. The protein localizes to lysosomal membranes and enables nicotinic acid adenine dinucleotide phosphate (NAADP) -induced calcium ion release from lysosome-related stores. This ubiquitously expressed gene has elevated expression in liver and kidney. Two common nonsynonymous SNPs in this gene strongly associate with blond versus brown hair pigmentation.[provided by RefSeq, Dec 2009]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <a href="mailto:techsupport@origene.com">techsupport@origene.com</a> . If you need a special design or shRNA sequence, please utilize our <a href="#">custom shRNA service</a> .



[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).