

## Product datasheet for **TL301025V**

### **TMEM126B Human shRNA Lentiviral Particle (Locus ID 55863)**

#### **Product data:**

<b>Product Type:</b>	shRNA Lentiviral Particles
<b>Product Name:</b>	TMEM126B Human shRNA Lentiviral Particle (Locus ID 55863)
<b>Locus ID:</b>	55863
<b>Synonyms:</b>	HT007; MC1DN29
<b>Vector:</b>	pGFP-C-shLenti (TR30023)
<b>Format:</b>	Lentiviral particles
<b>Components:</b>	TMEM126B - Human shRNA lentiviral particles (4 unique 29mer target-specific shRNA, 1 scramble control), 0.5 ml each, >10 <sup>7</sup> TU/ml.
<b>RefSeq:</b>	<a href="#">BC017574</a> , <a href="#">NM_001193537</a> , <a href="#">NM_001193538</a> , <a href="#">NM_001256546</a> , <a href="#">NM_001256547</a> , <a href="#">NM_018480</a> , <a href="#">NR_036465</a> , <a href="#">NM_001350393</a> , <a href="#">NM_001350394</a> , <a href="#">NM_001350395</a> , <a href="#">NM_001350396</a> , <a href="#">NR_146645</a> , <a href="#">NR_146646</a> , <a href="#">NM_018480.1</a> , <a href="#">NM_018480.2</a> , <a href="#">NM_018480.3</a> , <a href="#">NM_018480.4</a> , <a href="#">NM_001193538.1</a> , <a href="#">NM_001193538.2</a> , <a href="#">NM_001193537.1</a> , <a href="#">NM_001193537.2</a> , <a href="#">NM_001256547.1</a> , <a href="#">NM_001256546.1</a> , <a href="#">BC038933</a> , <a href="#">BC012065</a> , <a href="#">BC107900</a> , <a href="#">BM746060</a> , <a href="#">NM_001256547.2</a> , <a href="#">NM_001256546.2</a> , <a href="#">NM_001193537.3</a> , <a href="#">NM_001193538.3</a>
<b>UniProt ID:</b>	<a href="#">Q8IUX1</a>
<b>Summary:</b>	This gene encodes a mitochondrial transmembrane protein which is a component of the mitochondrial complex I assembly complex. The encoded protein serves as an assembly factor that is required for formation of the membrane arm of the complex. It interacts with NADH dehydrogenase [ubiquinone] 1 alpha subcomplex assembly factor 13. Naturally occurring mutations in this gene are associated with isolated complex I deficiency. A pseudogene of this gene has been defined on chromosome 9. [provided by RefSeq, Apr 2017]
<b>shRNA Design:</b>	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <a href="mailto:techsupport@origene.com">techsupport@origene.com</a> . If you need a special design or shRNA sequence, please utilize our <a href="#">custom shRNA service</a> .



[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).