

## Product datasheet for **TL301151**

### Testin (TES) Human shRNA Plasmid Kit (Locus ID 26136)

#### Product data:

Product Type:	shRNA Plasmids
Product Name:	Testin (TES) Human shRNA Plasmid Kit (Locus ID 26136)
Locus ID:	26136
Synonyms:	TESS; TESS-2
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	TES - Human, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 26136). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	<a href="#">NM_015641</a> , <a href="#">NM_152829</a> , <a href="#">NM_015641.1</a> , <a href="#">NM_015641.2</a> , <a href="#">NM_015641.3</a> , <a href="#">NM_152829.1</a> , <a href="#">NM_152829.2</a> , <a href="#">BC001451</a> , <a href="#">BC001451.1</a> , <a href="#">BC040208</a> , <a href="#">BC045580</a> , <a href="#">NM_152829.3</a> , <a href="#">NM_015641.4</a>
UniProt ID:	<a href="#">Q9UGI8</a>
Summary:	Cancer-associated chromosomal changes often involve regions containing fragile sites. This gene maps to a common fragile site on chromosome 7q31.2 designated FRA7G. This gene is similar to mouse Testin, a testosterone-responsive gene encoding a Sertoli cell secretory protein containing three LIM domains. LIM domains are double zinc-finger motifs that mediate protein-protein interactions between transcription factors, cytoskeletal proteins and signaling proteins. This protein is a negative regulator of cell growth and may act as a tumor suppressor. This scaffold protein may also play a role in cell adhesion, cell spreading and in the reorganization of the actin cytoskeleton. Multiple protein isoforms are encoded by transcript variants of this gene.[provided by RefSeq, Mar 2011]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <a href="mailto:techsupport@origene.com">techsupport@origene.com</a> . If you need a special design or shRNA sequence, please utilize our <a href="#">custom shRNA service</a> .


[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).