

## Product datasheet for **TL301196**

### **TBLR1 (TBL1XR1) Human shRNA Plasmid Kit (Locus ID 79718)**

#### **Product data:**

<b>Product Type:</b>	shRNA Plasmids
<b>Product Name:</b>	TBLR1 (TBL1XR1) Human shRNA Plasmid Kit (Locus ID 79718)
<b>Locus ID:</b>	79718
<b>Synonyms:</b>	C21; DC42; IRA1; MRD41; TBLR1
<b>Vector:</b>	pGFP-C-shLenti (TR30023)
<b>E. coli Selection:</b>	Chloramphenicol (34 ug/ml)
<b>Mammalian Cell Selection:</b>	Puromycin
<b>Format:</b>	Lentiviral plasmids
<b>Components:</b>	TBL1XR1 - Human, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 79718). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
<b>RefSeq:</b>	<a href="#">NM_001321193</a> , <a href="#">NM_001321194</a> , <a href="#">NM_001321195</a> , <a href="#">NM_024665</a> , <a href="#">NM_024665.1</a> , <a href="#">NM_024665.2</a> , <a href="#">NM_024665.3</a> , <a href="#">NM_024665.4</a> , <a href="#">NM_024665.5</a> , <a href="#">BC113421</a> , <a href="#">BC034712</a> , <a href="#">BC048298</a> , <a href="#">BC060320</a> , <a href="#">BC070134</a> , <a href="#">BC092400</a> , <a href="#">NM_024665.7</a>
<b>UniProt ID:</b>	<a href="#">Q9BZK7</a>
<b>Summary:</b>	This gene is a member of the WD40 repeat-containing gene family and shares sequence similarity with transducin (beta)-like 1X-linked (TBL1X). The protein encoded by this gene is thought to be a component of both nuclear receptor corepressor (N-CoR) and histone deacetylase 3 (HDAC 3) complexes, and is required for transcriptional activation by a variety of transcription factors. Mutations in these gene have been associated with some autism spectrum disorders, and one finding suggests that haploinsufficiency of this gene may be a cause of intellectual disability with dysmorphism. Mutations in this gene as well as recurrent translocations involving this gene have also been observed in some tumors. [provided by RefSeq, Mar 2016]
<b>shRNA Design:</b>	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <a href="mailto:techsupport@origene.com">techsupport@origene.com</a> . If you need a special design or shRNA sequence, please utilize our <a href="#">custom shRNA service</a> .



[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).