

## Product datasheet for **TL301248**

### Tapasin (TAPBP) Human shRNA Plasmid Kit (Locus ID 6892)

#### Product data:

Product Type:	shRNA Plasmids
Product Name:	Tapasin (TAPBP) Human shRNA Plasmid Kit (Locus ID 6892)
Locus ID:	6892
Synonyms:	NGS17; TAPA; TPN; TPSN
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	TAPBP - Human, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 6892). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	<a href="#">NM_003190</a> , <a href="#">NM_172208</a> , <a href="#">NM_172209</a> , <a href="#">NM_003190.1</a> , <a href="#">NM_003190.2</a> , <a href="#">NM_003190.3</a> , <a href="#">NM_003190.4</a> , <a href="#">NM_172208.1</a> , <a href="#">NM_172208.2</a> , <a href="#">NM_172209.1</a> , <a href="#">NM_172209.2</a> , <a href="#">BC064966</a> , <a href="#">BC080574</a> , <a href="#">NM_003190.5</a>
UniProt ID:	<a href="#">O15533</a>
Summary:	This gene encodes a transmembrane glycoprotein which mediates interaction between newly assembled major histocompatibility complex (MHC) class I molecules and the transporter associated with antigen processing (TAP), which is required for the transport of antigenic peptides across the endoplasmic reticulum membrane. This interaction is essential for optimal peptide loading on the MHC class I molecule. Up to four complexes of MHC class I and this protein may be bound to a single TAP molecule. This protein contains a C-terminal double-lysine motif (KKKAE) known to maintain membrane proteins in the endoplasmic reticulum. This gene lies within the major histocompatibility complex on chromosome 6. Alternative splicing results in three transcript variants encoding different isoforms. [provided by RefSeq, Jul 2008]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <a href="mailto:techsupport@origene.com">techsupport@origene.com</a> . If you need a special design or shRNA sequence, please utilize our <a href="#">custom shRNA service</a> .



[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).