

## Product datasheet for **TL301326V**

### Syntaxin 18 (STX18) Human shRNA Lentiviral Particle (Locus ID 53407)

#### Product data:

Product Type:	shRNA Lentiviral Particles
Product Name:	Syntaxin 18 (STX18) Human shRNA Lentiviral Particle (Locus ID 53407)
Locus ID:	53407
Synonyms:	Ufe1
Vector:	pGFP-C-shLenti (TR30023)
Format:	Lentiviral particles
Components:	STX18 - Human shRNA lentiviral particles (4 unique 29mer target-specific shRNA, 1 scramble control), 0.5 ml each, >10 <sup>7</sup> TU/ml.
RefSeq:	<a href="#">NM_016930</a> , <a href="#">NM_001346281</a> , <a href="#">NM_001346282</a> , <a href="#">NM_001346300</a> , <a href="#">NR_144414</a> , <a href="#">NM_016930.1</a> , <a href="#">NM_016930.2</a> , <a href="#">NM_016930.3</a> , <a href="#">BC014613</a> , <a href="#">BC014613.2</a> , <a href="#">BC071980</a> , <a href="#">NM_016930.4</a>
UniProt ID:	<a href="#">Q9P2W9</a>
Summary:	This gene encodes a member of the syntaxin family of soluble N-ethylmaleimide-sensitive factor attachment protein receptors (SNAREs) which is part of a membrane tethering complex that includes other SNAREs and several peripheral membrane proteins, and is involved in vesicular transport between the endoplasmic reticulum (ER) and the Golgi complex. The encoded protein is important for the organization of the smooth, rough, and exit site ER subdomains. A pseudogene of this gene has been identified. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Oct 2016]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <a href="mailto:techsupport@origene.com">techsupport@origene.com</a> . If you need a special design or shRNA sequence, please utilize our <a href="#">custom shRNA service</a> .



[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).