

## Product datasheet for **TL301357**

### Stabilin 2 (STAB2) Human shRNA Plasmid Kit (Locus ID 55576)

#### Product data:

Product Type:	shRNA Plasmids
Product Name:	Stabilin 2 (STAB2) Human shRNA Plasmid Kit (Locus ID 55576)
Locus ID:	55576
Synonyms:	FEEL2; FELE-2; FELL2; FEX2; HARE; SCARH1
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	STAB2 - Human, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 55576). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	<a href="#">NM_017564</a> , <a href="#">NM_017564.1</a> , <a href="#">NM_017564.2</a> , <a href="#">NM_017564.3</a> , <a href="#">NM_017564.4</a> , <a href="#">NM_017564.5</a> , <a href="#">NM_017564.6</a> , <a href="#">NM_017564.7</a> , <a href="#">NM_017564.8</a> , <a href="#">NM_017564.9</a> , <a href="#">BC172224</a>
UniProt ID:	<a href="#">Q8WWQ8</a>
Summary:	This gene encodes a large, transmembrane receptor protein which may function in angiogenesis, lymphocyte homing, cell adhesion, or receptor scavenging. The protein contains 7 fasciclin, 15 epidermal growth factor (EGF)-like, and 2 laminin-type EGF-like domains as well as a C-type lectin-like hyaluronan-binding Link module. The protein is primarily expressed on sinusoidal endothelial cells of liver, spleen, and lymph node. The receptor has been shown to bind and endocytose ligands such as hyaluronan, low density lipoprotein, Gram-positive and Gram-negative bacteria, and advanced glycosylation end products. Supporting its possible role as a scavenger receptor, the protein has been shown to cycle between the plasma membrane and lysosomes. [provided by RefSeq, Jul 2008]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <a href="mailto:techsupport@origene.com">techsupport@origene.com</a> . If you need a special design or shRNA sequence, please utilize our <a href="#">custom shRNA service</a> .



[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).