

Product datasheet for **TL301446**

SPAG11B Human shRNA Plasmid Kit (Locus ID 10407)

Product data:

Product Type:	shRNA Plasmids
Product Name:	SPAG11B Human shRNA Plasmid Kit (Locus ID 10407)
Locus ID:	10407
Synonyms:	EDDM2B; EP2; EP2C; EP2D; HE2; HE2C; SPAG11; SPAG11A
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	SPAG11B - Human, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 10407). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	NM_016512 , NM_058200 , NM_058201 , NM_058202 , NM_058203 , NM_058206 , NM_058207 , NM_016512.1 , NM_016512.2 , NM_016512.3 , NM_058202.1 , NM_058202.2 , NM_058200.1 , NM_058200.2 , NM_058207.1 , NM_058207.2 , NM_058203.1 , NM_058203.2 , NM_058201.1 , NM_058201.2 , NM_058206.1 , NM_058206.2 , NM_058206.3 , BC160085 , NM_058206.4 , NM_058207.3
UniProt ID:	Q6PDA7
Summary:	This gene encodes several androgen-dependent, epididymis-specific secretory proteins. The specific functions of these proteins have not been determined, but they are thought to be involved in sperm maturation. Some of the isoforms contain regions of similarity to beta-defensins, a family of antimicrobial peptides. The gene is located on chromosome 8p23 near the defensin gene cluster. Alternative splicing of this gene results in seven transcript variants encoding different isoforms. Two different N-terminal and five different C-terminal protein sequences are encoded by the splice variants. Two additional variants have been described, but their full length sequences have not been determined. [provided by RefSeq, Jul 2008]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com . If you need a special design or shRNA sequence, please utilize our custom shRNA service .



[View online »](#)

**Performance
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).