

Product datasheet for **TL301456**

SOX5 Human shRNA Plasmid Kit (Locus ID 6660)

Product data:

Product Type:	shRNA Plasmids
Product Name:	SOX5 Human shRNA Plasmid Kit (Locus ID 6660)
Locus ID:	6660
Synonyms:	L-SOX5; L-SOX5B; L-SOX5F; LAMSHF
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	SOX5 - Human, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 6660). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	BC014929 , NM_001261414 , NM_001261415 , NM_001330785 , NM_006940 , NM_152989 , NM_178010 , NM_178010.1 , NM_178010.2 , NM_178010.3 , NM_006940.1 , NM_006940.2 , NM_006940.3 , NM_006940.4 , NM_006940.5 , NM_152989.1 , NM_152989.2 , NM_152989.3 , NM_152989.4 , NM_001261414.1 , NM_001261414.2 , NM_001261415.1 , NM_001261415.2 , BC029220 , BC035227 , BC047665 , BC060773 , BC065532 , BM988193 , NM_152989.5 , NM_001261414.3 , NM_006940.6 , NM_178010.4
UniProt ID:	P35711
Summary:	This gene encodes a member of the SOX (SRY-related HMG-box) family of transcription factors involved in the regulation of embryonic development and in the determination of the cell fate. The encoded protein may act as a transcriptional regulator after forming a protein complex with other proteins. The encoded protein may play a role in chondrogenesis. A pseudogene of this gene is located on chromosome 8. Multiple transcript variants encoding distinct isoforms have been identified for this gene. [provided by RefSeq, Jul 2008]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com . If you need a special design or shRNA sequence, please utilize our custom shRNA service .



[View online »](#)

**Performance
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).