

## Product datasheet for **TL301474**

### SNURF Human shRNA Plasmid Kit (Locus ID 8926)

#### Product data:

Product Type:	shRNA Plasmids
Product Name:	SNURF Human shRNA Plasmid Kit (Locus ID 8926)
Locus ID:	8926
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	SNURF - Human, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 8926). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	<a href="#">NM_005678</a> , <a href="#">NM_022804</a> , <a href="#">NM_005678.2</a> , <a href="#">NM_005678.3</a> , <a href="#">NM_005678.4</a> , <a href="#">NM_022804.1</a> , <a href="#">NM_022804.2</a> , <a href="#">NM_022804.3</a> , <a href="#">NM_005678.1</a>
UniProt ID:	<a href="#">Q9Y675</a>
Summary:	This gene is located within the Prader-Willi Syndrome critical region on chromosome 15. Transcripts produced from this gene initiate at an imprinting center and are paternally-imprinted. These transcripts may be bicistronic and also encode SNRPN (small nuclear ribonucleoprotein polypeptide N) from a downstream open reading frame. The small protein represented by this gene is encoded by an evolutionarily-conserved upstream open reading frame and is localized to the nucleus. Extensive alternative splicing and promoter usage occurs in this region and the full-length nature of some of these transcripts has not been determined. Alterations in the imprinting center are associated with parental imprint switch failure, which may cause Angelman syndrome or Prader-Willi syndrome. [provided by RefSeq, Mar 2017]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <a href="mailto:techsupport@origene.com">techsupport@origene.com</a> . If you need a special design or shRNA sequence, please utilize our <a href="#">custom shRNA service</a> .



[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).