

## Product datasheet for **TL301502**

### SMC6L1 (SMC6) Human shRNA Plasmid Kit (Locus ID 79677)

#### Product data:

Product Type:	shRNA Plasmids
Product Name:	SMC6L1 (SMC6) Human shRNA Plasmid Kit (Locus ID 79677)
Locus ID:	79677
Synonyms:	hSMC6; SMC-6; SMC6L1
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	SMC6 - Human, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 79677). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	<a href="#">NM_001142286</a> , <a href="#">NM_024624</a> , <a href="#">NM_024624.1</a> , <a href="#">NM_024624.2</a> , <a href="#">NM_024624.3</a> , <a href="#">NM_024624.4</a> , <a href="#">NM_024624.5</a> , <a href="#">NM_001142286.1</a> , <a href="#">BC039828</a> , <a href="#">BC039828.1</a> , <a href="#">BC011389</a> , <a href="#">BC022998</a> , <a href="#">BC032675</a> , <a href="#">NM_001142286.2</a>
UniProt ID:	<a href="#">Q96SB8</a>
Summary:	Core component of the SMC5-SMC6 complex, a complex involved in DNA double-strand breaks by homologous recombination. The complex may promote sister chromatid homologous recombination by recruiting the SMC1-SMC3 cohesin complex to double-strand breaks. The complex is required for telomere maintenance via recombination in ALT (alternative lengthening of telomeres) cell lines and mediates sumoylation of shelterin complex (telosome) components which is proposed to lead to shelterin complex disassembly in ALT-associated PML bodies (APBs). Required for recruitment of telomeres to PML nuclear bodies. SMC5-SMC6 complex may prevent transcription of episomal DNA, such as circular viral DNA genome (PubMed:26983541).[UniProtKB/Swiss-Prot Function]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <a href="mailto:techsupport@origene.com">techsupport@origene.com</a> . If you need a special design or shRNA sequence, please utilize our <a href="#">custom shRNA service</a> .



[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).