

Product datasheet for **TL301620**

SLC26A9 Human shRNA Plasmid Kit (Locus ID 115019)

Product data:

Product Type:	shRNA Plasmids
Product Name:	SLC26A9 Human shRNA Plasmid Kit (Locus ID 115019)
Locus ID:	115019
Synonyms:	anion transporter/exchanger-9; OTTHUMP00000034236; OTTHUMP00000034326; solute carrier family 26, member 9
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	SLC26A9 - Human, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 115019). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	NM_001142600 , NM_052934 , NM_134325 , NM_052934.2 , NM_052934.3 , NM_134325.1 , NM_134325.2 , NM_001142600.1 , BC151208 , BC037978 , BC136538 , BC144233 , NM_052934.4
UniProt ID:	Q7LBE3
Summary:	This gene is one member of a family of sulfate/anion transporter genes. Family members are well conserved in their genomic (number and size of exons) and protein (aa length among species) structures yet have markedly different tissue expression patterns. The product of this gene is a highly selective chloride ion channel regulated by WNK kinases. Alternative splicing results in multiple transcript variants encoding differing isoforms.[provided by RefSeq, Dec 2008]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com . If you need a special design or shRNA sequence, please utilize our custom shRNA service .



[View online »](#)

**Performance
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).