

Product datasheet for **TL301682V**

LY108 (SLAMF6) Human shRNA Lentiviral Particle (Locus ID 114836)

Product data:

Product Type:	shRNA Lentiviral Particles
Product Name:	LY108 (SLAMF6) Human shRNA Lentiviral Particle (Locus ID 114836)
Locus ID:	114836
Synonyms:	CD352; KALI; KAL1b; Ly108; NTB-A; NTBA; SF2000
Vector:	pGFP-C-shLenti (TR30023)
Format:	Lentiviral particles
Components:	SLAMF6 - Human shRNA lentiviral particles (4 unique 29mer target-specific shRNA, 1 scramble control), 0.5 ml each, >10 ⁷ TU/ml.
RefSeq:	NM_001184714 , NM_001184715 , NM_001184716 , NM_052931 , NM_052931.1 , NM_052931.2 , NM_052931.4 , NM_001184716.1 , NM_001184715.1 , NM_001184714.1 , BC113893 , BC090928 , BC114495 , NM_052931.5 , NM_001184714.2 , NM_001184715.2 , NM_001184716.2
UniProt ID:	Q96DU3
Summary:	The protein encoded by this gene is a type I transmembrane protein, belonging to the CD2 subfamily of the immunoglobulin superfamily. This encoded protein is expressed on Natural killer (NK), T, and B lymphocytes. It undergoes tyrosine phosphorylation and associates with the Src homology 2 domain-containing protein (SH2D1A) as well as with SH2 domain-containing phosphatases (SHPs). It functions as a coreceptor in the process of NK cell activation. It can also mediate inhibitory signals in NK cells from X-linked lymphoproliferative patients. Alternative splicing results in multiple transcript variants encoding distinct isoforms. [provided by RefSeq, May 2010]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com . If you need a special design or shRNA sequence, please utilize our custom shRNA service .



[View online »](#)

**Performance
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).