

Product datasheet for **TL301926**

ROBO1 Human shRNA Plasmid Kit (Locus ID 6091)

Product data:

Product Type:	shRNA Plasmids
Product Name:	ROBO1 Human shRNA Plasmid Kit (Locus ID 6091)
Locus ID:	6091
Synonyms:	axon guidance receptor; DUTT1; FLJ21882; MGC131599; MGC133277; roundabout, axon guidance receptor, homolog 1 (Drosophila); roundabout 1; SAX3
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	ROBO1 - Human, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 6091). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	NM_001145844 , NM_001145845 , NM_002941 , NM_133631 , NM_133631.1 , NM_133631.2 , NM_133631.3 , NM_002941.1 , NM_002941.2 , NM_002941.3 , NM_001145845.1 , NM_001145844.1 , BC157861 , BC115022 , BC001969 , BC029676 , BC112336 , BC115020 , BC171855 , NM_002941.4
UniProt ID:	Q9Y6N7
Summary:	Bilateral symmetric nervous systems have special midline structures that establish a partition between the two mirror image halves. Some axons project toward and across the midline in response to long-range chemoattractants emanating from the midline. The product of this gene is a member of the immunoglobulin gene superfamily and encodes an integral membrane protein that functions in axon guidance and neuronal precursor cell migration. This receptor is activated by SLIT-family proteins, resulting in a repulsive effect on glioma cell guidance in the developing brain. A related gene is located at an adjacent region on chromosome 3. Multiple transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Mar 2009]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com . If you need a special design or shRNA sequence, please utilize our custom shRNA service .




[View online »](#)


**Performance
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.


For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).

Product images:


GFP signal was observed under microscope at 48 hours after transduction of TL301926A virus into HEK293 cells. TL301926A virus was prepared using lenti-shRNA TL301926A and [TR30037] packaging kit.



GFP signal was observed under microscope at 48 hours after transduction of TL301926B virus into HEK293 cells. TL301926B virus was prepared using lenti-shRNA TL301926B and [TR30037] packaging kit.



GFP signal was observed under microscope at 48 hours after transduction of [TL301926C] virus into HEK293 cells. [TL301926C] virus was prepared using lenti-shRNA [TL301926C] and [TR30037] packaging kit.



GFP signal was observed under microscope at 48 hours after transduction of [TL301926D] virus into HEK293 cells. [TL301926D] virus was prepared using lenti-shRNA [TL301926D] and [TR30037] packaging kit.