

## Product datasheet for **TL301962**

### **RNF146 Human shRNA Plasmid Kit (Locus ID 81847)**

#### **Product data:**

Product Type:	shRNA Plasmids
Product Name:	RNF146 Human shRNA Plasmid Kit (Locus ID 81847)
Locus ID:	81847
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	RNF146 - Human, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 81847). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	<a href="#">NM_001242844</a> , <a href="#">NM_001242845</a> , <a href="#">NM_001242846</a> , <a href="#">NM_001242847</a> , <a href="#">NM_001242848</a> , <a href="#">NM_001242849</a> , <a href="#">NM_001242850</a> , <a href="#">NM_001242851</a> , <a href="#">NM_001242852</a> , <a href="#">NM_030963</a> , <a href="#">NR_045102</a> , <a href="#">NR_045103</a> , <a href="#">NM_030963.1</a> , <a href="#">NM_030963.2</a> , <a href="#">NM_030963.3</a> , <a href="#">NM_001242844.1</a> , <a href="#">NM_001242845.1</a> , <a href="#">NM_001242846.1</a> , <a href="#">NM_001242852.1</a> , <a href="#">NM_001242849.1</a> , <a href="#">NM_001242850.1</a> , <a href="#">NM_001242851.1</a> , <a href="#">BC008235</a> , <a href="#">BC008235.1</a> , <a href="#">NM_001242850.2</a>
UniProt ID:	<a href="#">Q9NTX7</a>



[View online »](#)

<b>Summary:</b>	<p>E3 ubiquitin-protein ligase that specifically binds poly-ADP-ribosylated (PARsylated) proteins and mediates their ubiquitination and subsequent degradation. May regulate many important biological processes, such as cell survival and DNA damage response. Acts as an activator of the Wnt signaling pathway by mediating the ubiquitination of PARsylated AXIN1 and AXIN2, 2 key components of the beta-catenin destruction complex. Acts in cooperation with tankyrase proteins (TNKS and TNKS2), which mediate PARsylation of target proteins AXIN1, AXIN2, BLZF1, CASC3, TNKS and TNKS2. Recognizes and binds tankyrase-dependent PARsylated proteins via its WWE domain and mediates their ubiquitination, leading to their degradation. Different ubiquitin linkage types have been observed: TNKS2 undergoes ubiquitination at 'Lys-48' and 'Lys-63', while AXIN1 is only ubiquitinated at 'Lys-48'. May regulate TNKS and TNKS2 subcellular location, preventing aggregation at a centrosomal location. Neuroprotective protein. Protects the brain against N-methyl-D-aspartate (NMDA) receptor-mediated glutamate excitotoxicity and ischemia, by interfering with PAR-induced cell death, called parthanatos. Prevents nuclear translocation of AIFM1 in a PAR-binding dependent manner. Does not affect PARP1 activation (By similarity). Protects against cell death induced by DNA damaging agents, such as N-methyl-N-nitro-N-nitrosoguanidine (MNNG) and rescues cells from G1 arrest. Promotes cell survival after gamma-irradiation. Facilitates DNA repair.[UniProtKB/Swiss-Prot Function]</p>
<b>shRNA Design:</b>	<p>These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <a href="mailto:techsupport@origene.com">techsupport@origene.com</a>. If you need a special design or shRNA sequence, please utilize our <a href="#">custom shRNA service</a>.</p>
<b>Performance Guaranteed:</b>	<p>OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.</p> <p>For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at <a href="mailto:techsupport@origene.com">techsupport@origene.com</a>. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).</p>