

## Product datasheet for **TL302006**

### RHBG Human shRNA Plasmid Kit (Locus ID 57127)

#### Product data:

Product Type:	shRNA Plasmids
Product Name:	RHBG Human shRNA Plasmid Kit (Locus ID 57127)
Locus ID:	57127
Synonyms:	SLC42A2
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	RHBG - Human, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 57127). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	<a href="#">NM_001256395</a> , <a href="#">NM_001256396</a> , <a href="#">NM_020407</a> , <a href="#">NR_026549</a> , <a href="#">NR_046115</a> , <a href="#">NR_146763</a> , <a href="#">NR_146764</a> , <a href="#">NR_146765</a> , <a href="#">NM_020407.1</a> , <a href="#">NM_020407.2</a> , <a href="#">NM_020407.3</a> , <a href="#">NM_020407.4</a> , <a href="#">NM_001256395.1</a> , <a href="#">NM_001256396.1</a> , <a href="#">BC065218</a> , <a href="#">NM_001256395.2</a> , <a href="#">NM_001256396.2</a> , <a href="#">NM_020407.5</a>
UniProt ID:	<a href="#">Q9H310</a>
Summary:	This gene encodes one of two non-erythroid members of the Rhesus (Rh) protein family. Non-erythroid Rh protein family members are mainly expressed in the kidney and belong to the methylammonium-ammonium permease/ammonia transporters superfamily. All Rh family proteins are predicted to be transmembrane proteins with 12 membrane spanning domains and intracytoplasmic N- and C-termini. Alternative splicing of this gene results in multiple transcript variants. [provided by RefSeq, Jan 2012]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <a href="mailto:techsupport@origene.com">techsupport@origene.com</a> . If you need a special design or shRNA sequence, please utilize our <a href="#">custom shRNA service</a> .



[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).