

Product datasheet for **TL302024**

RFWD3 Human shRNA Plasmid Kit (Locus ID 55159)

Product data:


Product Type:	shRNA Plasmids
Product Name:	RFWD3 Human shRNA Plasmid Kit (Locus ID 55159)
Locus ID:	55159
Synonyms:	FANCW; RNF201
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	RFWD3 - Human, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 55159). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	NM_018124 , NM_018124.1 , NM_018124.2 , NM_018124.3 , BC059371 , BC059371.1 , BC002574 , BM467919 , NM_001370536 , NM_001370534 , NM_001370535 , NM_001370539 , NM_001370540 , NM_001370543 , NM_001370537 , NM_001370542
UniProt ID:	Q6PCD5




[View online »](#)

Summary:	<p>E3 ubiquitin-protein ligase required for the repair of DNA interstrand cross-links (ICL) in response to DNA damage (PubMed:21504906, PubMed:21558276, PubMed:26474068, PubMed:28575657, PubMed:28575658). Plays a key role in RPA-mediated DNA damage signaling and repair (PubMed:21504906, PubMed:21558276, PubMed:26474068, PubMed:28575657, PubMed:28575658, PubMed:28691929). Acts by mediating ubiquitination of the RPA complex (RPA1, RPA2 and RPA3 subunits) and RAD51 at stalled replication forks, leading to remove them from DNA damage sites and promote homologous recombination (PubMed:26474068, PubMed:28575657, PubMed:28575658). Also mediates the ubiquitination of p53/TP53 in the late response to DNA damage, and acts as a positive regulator of p53/TP53 stability, thereby regulating the G1/S DNA damage checkpoint (PubMed:20173098). May act by catalyzing the formation of short polyubiquitin chains on p53/TP53 that are not targeted to the proteasome (PubMed:20173098). In response to ionizing radiation, interacts with MDM2 and enhances p53/TP53 ubiquitination, possibly by restricting MDM2 from extending polyubiquitin chains on ubiquitinated p53/TP53 (PubMed:20173098). [UniProtKB/Swiss-Prot Function]</p>
shRNA Design:	<p>These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com. If you need a special design or shRNA sequence, please utilize our custom shRNA service.</p>
Performance Guaranteed:	<p>OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.</p> <p>For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).</p>


Product images:




GFP signal was observed under microscope at 48 hours after transduction of TL302024A virus into HEK293 cells. TL302024A virus was prepared using lenti-shRNA TL302024A and [TR30037] packaging kit.



GFP signal was observed under microscope at 48 hours after transduction of TL302024B virus into HEK293 cells. TL302024B virus was prepared using lenti-shRNA TL302024B and [TR30037] packaging kit.



GFP signal was observed under microscope at 48 hours after transduction of [TL302024C] virus into HEK293 cells. [TL302024C] virus was prepared using lenti-shRNA [TL302024C] and [TR30037] packaging kit.



GFP signal was observed under microscope at 48 hours after transduction of [TL302024D] virus into HEK293 cells. [TL302024D] virus was prepared using lenti-shRNA [TL302024D] and [TR30037] packaging kit.