

Product datasheet for **TL302087**

RBED1 (ELMOD3) Human shRNA Plasmid Kit (Locus ID 84173)

Product data:

Product Type:	shRNA Plasmids
Product Name:	RBED1 (ELMOD3) Human shRNA Plasmid Kit (Locus ID 84173)
Locus ID:	84173
Synonyms:	DFNB88; LST3; RBED1; RBM29
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	ELMOD3 - Human, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 84173). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	NM_001135021 , NM_001135022 , NM_001135023 , NM_001329791 , NM_001329792 , NM_001329793 , NM_032213 , NR_138131 , NR_138132 , NR_138133 , NM_032213.1 , NM_032213.2 , NM_032213.3 , NM_032213.4 , NM_001135021.1 , NM_001135022.1 , NM_001135023.1 , BC018666 , BC018666.2 , BC112324 , BC001942 , BC010991 , BC054504 , NM_001135021.2 , NM_001135023.2 , NM_001135022.2
UniProt ID:	Q96FG2
Summary:	This gene encodes a member of the engulfment and cell motility family of GTPase-activating proteins that regulate Arf GTPase proteins. Members of this family are defined by a conserved engulfment and cell motility domain. In rat cochlea, the encoded protein is found in stereocilia, kinocilia and cuticular plate of developing hair cells suggesting a function for this protein in cochlear sensory cells. An allelic variant of this family has been associated with autosomal recessive nonsyndromic deafness-88 in humans. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Jul 2016]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com . If you need a special design or shRNA sequence, please utilize our custom shRNA service .



[View online »](#)

**Performance
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).