

Product datasheet for **TL302117**

RANBP17 Human shRNA Plasmid Kit (Locus ID 64901)

Product data:

Product Type:	shRNA Plasmids
Product Name:	RANBP17 Human shRNA Plasmid Kit (Locus ID 64901)
Locus ID:	64901
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	RANBP17 - Human, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 64901). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	NM_022897 , NM_022897.1 , NM_022897.2 , NM_022897.3 , NM_022897.4 , BC030953 , BC148559 , NM_022897.5
UniProt ID:	Q9H2T7
Summary:	The transport of protein and large RNAs through the nuclear pore complexes (NPC) is an energy-dependent and regulated process. The import of proteins with a nuclear localization signal (NLS) is accomplished by recognition of one or more clusters of basic amino acids by the importin-alpha/beta complex; see MIM 600685 and MIM 602738. The small GTPase RAN (MIM 601179) plays a key role in NLS-dependent protein import. RAN-binding protein-17 is a member of the importin-beta superfamily of nuclear transport receptors.[supplied by OMIM, Jul 2002]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com . If you need a special design or shRNA sequence, please utilize our custom shRNA service .




[View online »](#)


**Performance
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.


For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).

Product images:


GFP signal was observed under microscope at 48 hours after transduction of TL302117A virus into HEK293 cells. TL302117A virus was prepared using lenti-shRNA TL302117A and [TR30037] packaging kit.



GFP signal was observed under microscope at 48 hours after transduction of TL302117B virus into HEK293 cells. TL302117B virus was prepared using lenti-shRNA TL302117B and [TR30037] packaging kit.



GFP signal was observed under microscope at 48 hours after transduction of [TL302117C] virus into HEK293 cells. [TL302117C] virus was prepared using lenti-shRNA [TL302117C] and [TR30037] packaging kit.



GFP signal was observed under microscope at 48 hours after transduction of [TL302117D] virus into HEK293 cells. [TL302117D] virus was prepared using lenti-shRNA [TL302117D] and [TR30037] packaging kit.