

## Product datasheet for **TL302140**

### **RAB6A Human shRNA Plasmid Kit (Locus ID 5870)**

#### **Product data:**

Product Type:	shRNA Plasmids
Product Name:	RAB6A Human shRNA Plasmid Kit (Locus ID 5870)
Locus ID:	5870
Synonyms:	RAB6
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	RAB6A - Human, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 5870). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	<a href="#">NM_001243718</a> , <a href="#">NM_001243719</a> , <a href="#">NM_002869</a> , <a href="#">NM_198896</a> , <a href="#">NM_002869.1</a> , <a href="#">NM_002869.2</a> , <a href="#">NM_002869.3</a> , <a href="#">NM_002869.4</a> , <a href="#">NM_198896.1</a> , <a href="#">NM_001243718.1</a> , <a href="#">NM_001243719.1</a> , <a href="#">BC003617</a> , <a href="#">BC003617.2</a> , <a href="#">BC096818</a> , <a href="#">BC096818.1</a> , <a href="#">BC034996</a> , <a href="#">BC068486</a> , <a href="#">NM_002869.5</a> , <a href="#">NM_198896.2</a> , <a href="#">NM_001243718.2</a>
UniProt ID:	<a href="#">P20340</a>
Summary:	This gene encodes a member of the RAB family, which belongs to the small GTPase superfamily. GTPases of the RAB family bind to various effectors to regulate the targeting and fusion of transport carriers to acceptor compartments. This protein is located at the Golgi apparatus, which regulates trafficking in both a retrograde (from early endosomes and Golgi to the endoplasmic reticulum) and an anterograde (from the Golgi to the plasma membrane) directions. Myosin II is an effector of this protein in these processes. This protein is also involved in assembly of human cytomegalovirus (HCMV) by interacting with the cellular protein Bicaudal D1, which interacts with the HCMV virion tegument protein, pp150. Multiple alternatively spliced transcript variants encoding different isoforms have been identified. [provided by RefSeq, Aug 2011]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <a href="mailto:techsupport@origene.com">techsupport@origene.com</a> . If you need a special design or shRNA sequence, please utilize our <a href="#">custom shRNA service</a> .



[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).