

Product datasheet for **TL302321V**

PPP1R9B Human shRNA Lentiviral Particle (Locus ID 84687)

Product data:

Product Type:	shRNA Lentiviral Particles
Product Name:	PPP1R9B Human shRNA Lentiviral Particle (Locus ID 84687)
Locus ID:	84687
Synonyms:	PPP1R6; PPP1R9; SPINO; Spn
Vector:	pGFP-C-shLenti (TR30023)
Format:	Lentiviral particles
Components:	PPP1R9B - Human shRNA lentiviral particles (4 unique 29mer target-specific shRNA, 1 scramble control), 0.5 ml each, >10 ⁷ TU/ml.
RefSeq:	NM_032595 , NM_032595.1 , NM_032595.2 , NM_032595.3 , NM_032595.4 , BC008564 , BC011232 , BC016162 , BM991997 , NM_032595.5
UniProt ID:	Q96SB3
Summary:	This gene encodes a scaffold protein that functions as a regulatory subunit of protein phosphatase 1a. Expression of this gene is particularly high in dendritic spines, suggesting that the encoded protein may play a role in receiving signals from the central nervous system. The encoded protein has putative tumor suppressor function and decreased expression has been observed in tumors. [provided by RefSeq, Feb 2014]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com . If you need a special design or shRNA sequence, please utilize our custom shRNA service .



[View online »](#)

**Performance
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).