

## Product datasheet for **TL302339V**

### PPM1B Human shRNA Lentiviral Particle (Locus ID 5495)

#### Product data:

Product Type:	shRNA Lentiviral Particles
Product Name:	PPM1B Human shRNA Lentiviral Particle (Locus ID 5495)
Locus ID:	5495
Synonyms:	PP2C-beta-X; PP2CB; PP2CBETA; PPC2BETAX
Vector:	pGFP-C-shLenti (TR30023)
Format:	Lentiviral particles
Components:	PPM1B - Human shRNA lentiviral particles (4 unique 29mer target-specific shRNA, 1 scramble control), 0.5 ml each, >10 <sup>7</sup> TU/ml.
RefSeq:	<a href="#">NM_001033556</a> , <a href="#">NM_001033557</a> , <a href="#">NM_002706</a> , <a href="#">NM_177968</a> , <a href="#">NM_177969</a> , <a href="#">NM_177969.1</a> , <a href="#">NM_177969.2</a> , <a href="#">NM_177969.3</a> , <a href="#">NM_002706.1</a> , <a href="#">NM_002706.2</a> , <a href="#">NM_002706.3</a> , <a href="#">NM_002706.4</a> , <a href="#">NM_002706.5</a> , <a href="#">NM_001033556.1</a> , <a href="#">NM_001033557.1</a> , <a href="#">NM_001033557.2</a> , <a href="#">NM_177968.1</a> , <a href="#">NM_177968.2</a> , <a href="#">NM_177968.3</a> , <a href="#">BC064381</a> , <a href="#">BC064381.1</a> , <a href="#">BC012002</a> , <a href="#">NM_002706.6</a> , <a href="#">NM_177968.4</a> , <a href="#">NM_001033557.3</a>
UniProt ID:	<a href="#">O75688</a>
Summary:	The protein encoded by this gene is a member of the PP2C family of Ser/Thr protein phosphatases. PP2C family members are known to be negative regulators of cell stress response pathways. This phosphatase has been shown to dephosphorylate cyclin-dependent kinases (CDKs), and thus may be involved in cell cycle control. Overexpression of this phosphatase is reported to cause cell-growth arrest or cell death. Alternative splicing results in multiple transcript variants encoding different isoforms. Additional transcript variants have been described, but currently do not represent full-length sequences. [provided by RefSeq, Jul 2008]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <a href="mailto:techsupport@origene.com">techsupport@origene.com</a> . If you need a special design or shRNA sequence, please utilize our <a href="#">custom shRNA service</a> .



[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).