

Product datasheet for **TL302394**

PNPLA1 Human shRNA Plasmid Kit (Locus ID 285848)

Product data:

Product Type:	shRNA Plasmids
Product Name:	PNPLA1 Human shRNA Plasmid Kit (Locus ID 285848)
Locus ID:	285848
Synonyms:	dj50j22.1; FLJ38755; FLJ38755, MGC119319, MGC119320, MGC119321, dj50j22.1; MGC119319; MGC119320; MGC119321; OTTHUMP00000016289; OTTHUMP00000039681; patatin-like phospholipase domain containing 1
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	PNPLA1 - Human, 4 unique 29mer shRNA constructs in lentiviral GFP vector (Gene ID = 285848). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	NM_001039725 , NM_001145716 , NM_001145717 , NM_173676 , NM_173676.1 , NM_173676.2 , NM_001145716.1 , NM_001145716.2 , NM_001145717.1 , NM_001039725.1 , BC103906 , BC103905 , BC103907
UniProt ID:	Q8N8W4
Summary:	The protein encoded by this gene belongs to the patatin-like phospholipase (PNPLA) family, which is characterized by the presence of a highly conserved patatin domain. PNPLA family members have diverse lipolytic and acyltransferase activities, and are key elements in lipid metabolism. While other members of this family have been well characterized, the function of this gene remained an enigma. However, recent studies show that this gene is expressed in the skin epidermal keratinocytes, and has a role in glycerophospholipid metabolism in the cutaneous barrier. Consistent with these observations, mutations in this gene are associated with ichthyosis in human (autosomal recessive congenital ichthyoses, ARCI) and dog. Alternatively spliced transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Jun 2012]


[View online »](#)

shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com . If you need a special design or shRNA sequence, please utilize our custom shRNA service .
Performance Guaranteed:	OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples. For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).