

Product datasheet for **TL302404V**

Plexin A3 (PLXNA3) Human shRNA Lentiviral Particle (Locus ID 55558)

Product data:

Product Type:	shRNA Lentiviral Particles
Product Name:	Plexin A3 (PLXNA3) Human shRNA Lentiviral Particle (Locus ID 55558)
Locus ID:	55558
Synonyms:	6.3; HSSEXGENE; PLXN3; PLXN4; XAP-6
Vector:	pGFP-C-shLenti (TR30023)
Format:	Lentiviral particles
Components:	PLXNA3 - Human shRNA lentiviral particles (4 unique 29mer target-specific shRNA, 1 scramble control), 0.5 ml each, >10 ⁷ TU/ml.
RefSeq:	NM_017514 , NM_017514.1 , NM_017514.2 , NM_017514.3 , NM_017514.4 , BC156373 , BC172494
UniProt ID:	P51805
Summary:	This gene encodes a member of the plexin class of proteins. The encoded protein is a class 3 semaphorin receptor, and may be involved in cytoskeletal remodeling and as well as apoptosis. Studies of a similar gene in zebrafish suggest that it is important for axon pathfinding in the developing nervous system. This gene may be associated with tumor progression. [provided by RefSeq, Aug 2013]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com . If you need a special design or shRNA sequence, please utilize our custom shRNA service .



[View online »](#)

**Performance
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).