

Product datasheet for **TL302493V**

PIAS1 Human shRNA Lentiviral Particle (Locus ID 8554)

Product data:

Product Type:	shRNA Lentiviral Particles
Product Name:	PIAS1 Human shRNA Lentiviral Particle (Locus ID 8554)
Locus ID:	8554
Synonyms:	DDXBP1; GBP; GU/RH-II; ZMIZ3
Vector:	pGFP-C-shLenti (TR30023)
Format:	Lentiviral particles
Components:	PIAS1 - Human shRNA lentiviral particles (4 unique 29mer target-specific shRNA, 1 scramble control), 0.5 ml each, >10 ⁷ TU/ml.
RefSeq:	NM_001320687 , NM_016166 , NM_016166.1 , NM_016166.2 , BC053685 , BC104646 , BC118587 , BC121797 , NM_016166.3
UniProt ID:	O75925
Summary:	This gene encodes a member of the protein inhibitor of activated STAT (PIAS) family. PIAS proteins function as SUMO E3 ligases and play important roles in many cellular processes by mediating the sumoylation of target proteins. This protein plays a central role as a transcriptional coregulator of numerous cellular pathways including the STAT1 and nuclear factor kappaB pathways. Alternate splicing results in multiple transcript variants. [provided by RefSeq, Mar 2016]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com . If you need a special design or shRNA sequence, please utilize our custom shRNA service .




[View online »](#)


**Performance
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).

Product images:

GFP signal was observed under microscope at 48 hours after transduction of [TL302493C] virus into HEK293 cells. [TL302493C] virus was prepared using lenti-shRNA [TL302493C] and [TR30037] packaging kit.



GFP signal was observed under microscope at 48 hours after transduction of [TL302493D] virus into HEK293 cells. [TL302493D] virus was prepared using lenti-shRNA [TL302493D] and [TR30037] packaging kit.