

## Product datasheet for **TL302510**

### PHF17 (JADE1) Human shRNA Plasmid Kit (Locus ID 79960)

#### Product data:

|                           |  |
|---------------------------|--|
| Product Type:             | shRNA Plasmids   |
| Product Name:             | PHF17 (JADE1) Human shRNA Plasmid Kit (Locus ID 79960)   |
| Locus ID:                 | 79960  |
| Synonyms:                 | PHF17  |
| Vector:                   | pGFP-C-shLenti (TR30023)   |
| E. coli Selection:        | Chloramphenicol (34 ug/ml)   |
| Mammalian Cell Selection: | Puromycin  |
| Format:                   | Lentiviral plasmids  |
| Components:               | JADE1 - Human, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 79960).<br>5µg purified plasmid DNA per construct<br>29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.   |
| RefSeq:                   | <a href="#">NM_001287437</a> , <a href="#">NM_001287439</a> , <a href="#">NM_001287440</a> , <a href="#">NM_001287441</a> , <a href="#">NM_001287442</a> ,<br><a href="#">NM_001287443</a> , <a href="#">NM_024900</a> , <a href="#">NM_199320</a> , <a href="#">NM_024900.1</a> , <a href="#">NM_024900.2</a> , <a href="#">NM_024900.3</a> ,<br><a href="#">NM_024900.4</a> , <a href="#">NM_199320.2</a> , <a href="#">NM_199320.3</a> , <a href="#">NM_001287441.1</a> , <a href="#">NM_001287437.1</a> ,<br><a href="#">NM_001287439.1</a> , <a href="#">NM_001287440.1</a> , <a href="#">NM_001287442.1</a> , <a href="#">NM_001287443.1</a> , <a href="#">BC032376</a> ,<br><a href="#">BC032376.1</a> , <a href="#">NM_199320.4</a> , <a href="#">NM_001287441.2</a> , <a href="#">NM_001287440.2</a> , <a href="#">NM_001287437.2</a> ,<br><a href="#">NM_001287439.2</a> , <a href="#">NM_024900.5</a> |
| UniProt ID:               | <a href="#">Q6IE81</a>   |
| Summary:                  | Component of the HBO1 complex which has a histone H4-specific acetyltransferase activity, a reduced activity toward histone H3 and is responsible for the bulk of histone H4 acetylation in vivo. Transcriptional coactivator, it may also promote acetylation of nucleosomal histone H4 by KAT5. Promotes apoptosis. May act as a renal tumor suppressor. Negatively regulates canonical Wnt signaling; at least in part, cooperates with NPHP4 in this function.<br>[UniProtKB/Swiss-Prot Function]  |
| shRNA Design:             | These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <a href="mailto:techsupport@origene.com">techsupport@origene.com</a> . If you need a special design or shRNA sequence, please utilize our <a href="#">custom shRNA service</a> .   |




[View online »](#)


**Performance  
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.


For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).

**Product images:**


GFP signal was observed under microscope at 48 hours after transduction of TL302510A virus into HEK293 cells. TL302510A virus was prepared using lenti-shRNA TL302510A and [TR30037] packaging kit.



GFP signal was observed under microscope at 48 hours after transduction of TL302510B virus into HEK293 cells. TL302510B virus was prepared using lenti-shRNA TL302510B and [TR30037] packaging kit.



GFP signal was observed under microscope at 48 hours after transduction of [TL302510C] virus into HEK293 cells. [TL302510C] virus was prepared using lenti-shRNA [TL302510C] and [TR30037] packaging kit.



GFP signal was observed under microscope at 48 hours after transduction of [TL302510D] virus into HEK293 cells. [TL302510D] virus was prepared using lenti-shRNA [TL302510D] and [TR30037] packaging kit.