

## Product datasheet for **TL302569**

### PDE11A Human shRNA Plasmid Kit (Locus ID 50940)

#### Product data:

Product Type:	shRNA Plasmids
Product Name:	PDE11A Human shRNA Plasmid Kit (Locus ID 50940)
Locus ID:	50940
Synonyms:	PPNAD2
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	PDE11A - Human, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 50940). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	<a href="#">NM_001077196</a> , <a href="#">NM_001077197</a> , <a href="#">NM_001077358</a> , <a href="#">NM_016953</a> , <a href="#">NM_001077358.1</a> , <a href="#">NM_016953.1</a> , <a href="#">NM_016953.2</a> , <a href="#">NM_016953.3</a> , <a href="#">NM_001077196.1</a> , <a href="#">NM_001077197.1</a> , <a href="#">BC114431</a> , <a href="#">BC112393</a> , <a href="#">NM_001077358.2</a> , <a href="#">NM_001077197.2</a> , <a href="#">NM_001077196.2</a> , <a href="#">NM_016953.4</a>
UniProt ID:	<a href="#">Q9HCR9</a>
Summary:	The 3',5'-cyclic nucleotides cAMP and cGMP function as second messengers in a wide variety of signal transduction pathways. 3',5'-cyclic nucleotide phosphodiesterases (PDEs) catalyze the hydrolysis of cAMP and cGMP to the corresponding 5'-monophosphates and provide a mechanism to downregulate cAMP and cGMP signaling. This gene encodes a member of the PDE protein superfamily. Mutations in this gene are a cause of Cushing disease and adrenocortical hyperplasia. Multiple transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Jul 2008]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <a href="mailto:techsupport@origene.com">techsupport@origene.com</a> . If you need a special design or shRNA sequence, please utilize our <a href="#">custom shRNA service</a> .



[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).