

## Product datasheet for **TL302639**

### PCDHA11 Human shRNA Plasmid Kit (Locus ID 56138)

#### Product data:

Product Type:	shRNA Plasmids
Product Name:	PCDHA11 Human shRNA Plasmid Kit (Locus ID 56138)
Locus ID:	56138
Synonyms:	CNR7; CNRN7; CNRS7; CRNR7; PCDH-ALPHA11
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	PCDHA11 - Human, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 56138). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	<a href="#">NM_018902</a> , <a href="#">NM_031861</a> , <a href="#">NM_031861.1</a> , <a href="#">NM_031861.2</a> , <a href="#">NM_018902.1</a> , <a href="#">NM_018902.2</a> , <a href="#">NM_018902.3</a> , <a href="#">NM_018902.4</a> , <a href="#">BC136706</a> , <a href="#">NM_018902.5</a>
UniProt ID:	<a href="#">Q9Y5I1</a>
Summary:	This gene is a member of the protocadherin alpha gene cluster, one of three related gene clusters tandemly linked on chromosome five that demonstrate an unusual genomic organization similar to that of B-cell and T-cell receptor gene clusters. The alpha gene cluster is composed of 15 cadherin superfamily genes related to the mouse CNR genes and consists of 13 highly similar and 2 more distantly related coding sequences. The tandem array of 15 N-terminal exons, or variable exons, are followed by downstream C-terminal exons, or constant exons, which are shared by all genes in the cluster. The large, uninterrupted N-terminal exons each encode six cadherin ectodomains while the C-terminal exons encode the cytoplasmic domain. These neural cadherin-like cell adhesion proteins are integral plasma membrane proteins that most likely play a critical role in the establishment and function of specific cell-cell connections in the brain. Alternative splicing has been observed and additional variants have been suggested but their full-length nature has yet to be determined. [provided by RefSeq, Jul 2008]




[View online »](#)

**shRNA Design:** These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact [techsupport@origene.com](mailto:techsupport@origene.com). If you need a special design or shRNA sequence, please utilize our [custom shRNA service](#).


**Performance Guaranteed:** OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).


### Product images:




GFP signal was observed under microscope at 48 hours after transduction of TL302639A virus into HEK293 cells. TL302639A virus was prepared using lenti-shRNA TL302639A and [TR30037] packaging kit.



GFP signal was observed under microscope at 48 hours after transduction of TL302639B virus into HEK293 cells. TL302639B virus was prepared using lenti-shRNA TL302639B and [TR30037] packaging kit.



GFP signal was observed under microscope at 48 hours after transduction of [TL302639C] virus into HEK293 cells. [TL302639C] virus was prepared using lenti-shRNA [TL302639C] and [TR30037] packaging kit.



GFP signal was observed under microscope at 48 hours after transduction of [TL302639D] virus into HEK293 cells. [TL302639D] virus was prepared using lenti-shRNA [TL302639D] and [TR30037] packaging kit.