

Product datasheet for **TL302664**

PAX3 Human shRNA Plasmid Kit (Locus ID 5077)

Product data:

Product Type:	shRNA Plasmids
Product Name:	PAX3 Human shRNA Plasmid Kit (Locus ID 5077)
Locus ID:	5077
Synonyms:	CDHS; HUP2; WS1; WS3
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	PAX3 - Human, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 5077). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	NM_000438 , NM_001127366 , NM_013942 , NM_181457 , NM_181458 , NM_181459 , NM_181460 , NM_181461 , NM_013942.1 , NM_013942.2 , NM_013942.3 , NM_013942.4 , NM_181457.1 , NM_181457.2 , NM_181457.3 , NM_181460.1 , NM_181460.2 , NM_181460.3 , NM_181458.1 , NM_181458.2 , NM_181458.3 , NM_181459.1 , NM_181459.2 , NM_181459.3 , NM_181461.1 , NM_181461.2 , NM_181461.3 , NM_000438.1 , NM_000438.2 , NM_000438.3 , NM_000438.4 , NM_000438.5 , NM_001127366.1 , NM_001127366.2 , BC101299 , BC101301 , BC063547 , BC101300 , BC101302 , BC114363 , NM_181459.4 , NM_000438.6 , NM_013942.5 , NM_181461.4 , NM_001127366.3 , NM_181460.4 , NM_181458.4
UniProt ID:	P23760
Summary:	This gene is a member of the paired box (PAX) family of transcription factors. Members of the PAX family typically contain a paired box domain and a paired-type homeodomain. These genes play critical roles during fetal development. Mutations in paired box gene 3 are associated with Waardenburg syndrome, craniofacial-deafness-hand syndrome, and alveolar rhabdomyosarcoma. The translocation t(2;13)(q35;q14), which represents a fusion between PAX3 and the forkhead gene, is a frequent finding in alveolar rhabdomyosarcoma. Alternative splicing results in transcripts encoding isoforms with different C-termini. [provided by RefSeq, Jul 2008]



[View online »](#)

- shRNA Design:** These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com. If you need a special design or shRNA sequence, please utilize our [custom shRNA service](#).
- Performance Guaranteed:** OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.
- For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).