

Product datasheet for **TL302698**

PANK1 Human shRNA Plasmid Kit (Locus ID 53354)

Product data:

Product Type:	shRNA Plasmids
Product Name:	PANK1 Human shRNA Plasmid Kit (Locus ID 53354)
Locus ID:	53354
Synonyms:	PANK
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	PANK1 - Human, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 53354). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	<u>NM_138316</u> , <u>NM_148977</u> , <u>NM_148978</u> , <u>NM_148978.1</u> , <u>NM_148978.2</u> , <u>NM_138316.1</u> , <u>NM_138316.2</u> , <u>NM_138316.3</u> , <u>NM_148977.1</u> , <u>NM_148977.2</u> , <u>BC026296</u> , <u>BC036934</u> , <u>BC156381</u> , <u>BC172476</u> , <u>NM_138316.4</u> , <u>NM_148978.3</u>
UniProt ID:	<u>Q8TE04</u>
Summary:	This gene encodes a member of the pantothenate kinase family. Pantothenate kinases are key regulatory enzymes in the biosynthesis of coenzyme A (CoA). The encoded protein catalyzes the first and rate-limiting enzymatic reaction in CoA biosynthesis and is regulated by CoA through feedback inhibition. Alternatively spliced transcript variants encoding multiple isoforms have been observed for this gene. This gene and an intronic miRNA on the same strand are co-regulated by the tumor suppressor p53 (see PMID 20833636). [provided by RefSeq, Apr 2011]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com . If you need a special design or shRNA sequence, please utilize our custom shRNA service .



[View online »](#)

**Performance
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).