

Product datasheet for **TL302759V**

ORC4L (ORC4) Human shRNA Lentiviral Particle (Locus ID 5000)

Product data:

Product Type:	shRNA Lentiviral Particles
Product Name:	ORC4L (ORC4) Human shRNA Lentiviral Particle (Locus ID 5000)
Locus ID:	5000
Synonyms:	ORC4L; ORC4P
Vector:	pGFP-C-shLenti (TR30023)
Format:	Lentiviral particles
Components:	ORC4 - Human shRNA lentiviral particles (4 unique 29mer target-specific shRNA, 1 scramble control), 0.5 ml each, >10 ⁷ TU/ml.
RefSeq:	NM_001190879 , NM_001190881 , NM_001190882 , NM_002552 , NM_181741 , NM_181742 , NM_002552.1 , NM_002552.2 , NM_002552.3 , NM_002552.4 , NM_181742.1 , NM_181742.2 , NM_181742.3 , NM_181741.1 , NM_181741.2 , NM_181741.3 , NM_001190881.1 , NM_001190881.2 , NM_001190882.1 , NM_001190882.2 , NM_001190879.1 , NM_001190879.2 , BC014847 , BC014847.1 , BC005388 , BC016177 , BM468640 , BM926755 , NM_181741.4 , NM_002552.5
UniProt ID:	O43929
Summary:	The origin recognition complex (ORC) is a highly conserved six subunit protein complex essential for the initiation of the DNA replication in eukaryotic cells. Studies in yeast demonstrated that ORC binds specifically to origins of replication and serves as a platform for the assembly of additional initiation factors such as Cdc6 and Mcm proteins. This gene encodes a subunit of the ORC complex. Several alternatively spliced transcript variants, some of which encode the same protein, have been reported for this gene. [provided by RefSeq, Oct 2010]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com . If you need a special design or shRNA sequence, please utilize our custom shRNA service .


[View online »](#)

**Performance
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).