

## Product datasheet for **TL302824**

### OBSCN Human shRNA Plasmid Kit (Locus ID 84033)

#### Product data:

Product Type:	shRNA Plasmids
Product Name:	OBSCN Human shRNA Plasmid Kit (Locus ID 84033)
Locus ID:	84033
Synonyms:	ARHGEF30; UNC89
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	OBSCN - Human, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 84033). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	<a href="#">NM_001098623</a> , <a href="#">NM_001271223</a> , <a href="#">NM_052843</a> , <a href="#">NM_052843.1</a> , <a href="#">NM_052843.2</a> , <a href="#">NM_052843.3</a> , <a href="#">NM_001098623.1</a> , <a href="#">BC034746</a> , <a href="#">BC043916</a> , <a href="#">BC098118</a> , <a href="#">BC098262</a> , <a href="#">BC098292</a> , <a href="#">BC099731</a> , <a href="#">BC114382</a> , <a href="#">BC114461</a>
UniProt ID:	<a href="#">Q5VST9</a>
Summary:	The obscurin gene spans more than 150 kb, contains over 80 exons and encodes a protein of approximately 720 kDa. The encoded protein contains 68 Ig domains, 2 fibronectin domains, 1 calcium/calmodulin-binding domain, 1 RhoGEF domain with an associated PH domain, and 2 serine-threonine kinase domains. This protein belongs to the family of giant sacromeric signaling proteins that includes titin and nebulin, and may have a role in the organization of myofibrils during assembly and may mediate interactions between the sarcoplasmic reticulum and myofibrils. Alternatively spliced transcript variants encoding different isoforms have been identified. [provided by RefSeq, Jul 2008]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <a href="mailto:techsupport@origene.com">techsupport@origene.com</a> . If you need a special design or shRNA sequence, please utilize our <a href="#">custom shRNA service</a> .



[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).