

Product datasheet for **TL302852**

SDOS (NUDT16L1) Human shRNA Plasmid Kit (Locus ID 84309)

Product data:

Product Type:	shRNA Plasmids
Product Name:	SDOS (NUDT16L1) Human shRNA Plasmid Kit (Locus ID 84309)
Locus ID:	84309
Synonyms:	SDOS; TIRR
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	NUDT16L1 - Human, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 84309). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	NM_001193452 , NM_032349 , NM_032349.1 , NM_032349.2 , NM_032349.3 , NM_001193452.1 , BC006223 , BC006223.2 , NM_001370587 , NM_001370585 , NM_001370586 , NM_001370588
UniProt ID:	Q9BRJ7
Summary:	Key regulator of TP53BP1 required to stabilize TP53BP1 and regulate its recruitment to chromatin (PubMed:28241136). In absence of DNA damage, interacts with the tandem Tudor-like domain of TP53BP1, masking the region that binds histone H4 dimethylated at 'Lys-20' (H4K20me2), thereby preventing TP53BP1 recruitment to chromatin and maintaining TP53BP1 localization to the nucleus (PubMed:28241136). Following DNA damage, ATM-induced phosphorylation of TP53BP1 and subsequent recruitment of RIF1 leads to dissociate NUDT16L1/TIRR from TP53BP1, unmasking the tandem Tudor-like domain and allowing recruitment of TP53BP1 to DNA double strand breaks (DSBs) (PubMed:28241136). Binds U8 snoRNA (PubMed:18820299).[UniProtKB/Swiss-Prot Function]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com . If you need a special design or shRNA sequence, please utilize our custom shRNA service .



[View online »](#)

**Performance
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).