

## Product datasheet for **TL303018**

### NBPF3 Human shRNA Plasmid Kit (Locus ID 84224)

#### Product data:

Product Type:	shRNA Plasmids
Product Name:	NBPF3 Human shRNA Plasmid Kit (Locus ID 84224)
Locus ID:	84224
Synonyms:	AE2
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	NBPF3 - Human, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 84224). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	<a href="#">NM_001256416</a> , <a href="#">NM_001256417</a> , <a href="#">NM_001330381</a> , <a href="#">NM_032264</a> , <a href="#">NR_046176</a> , <a href="#">NM_032264.1</a> , <a href="#">NM_032264.2</a> , <a href="#">NM_032264.3</a> , <a href="#">NM_032264.4</a> , <a href="#">NM_001256417.2</a> , <a href="#">NM_001256416.1</a> , <a href="#">NM_001256416.2</a> , <a href="#">BC024011</a> , <a href="#">BC146688</a> , <a href="#">NM_001256416.4</a> , <a href="#">NM_001256417.4</a> , <a href="#">NM_032264.6</a>
UniProt ID:	<a href="#">Q9H094</a>
Summary:	This gene is a member of the neuroblastoma breakpoint family (NBPF) which consists of dozens of recently duplicated genes primarily located in segmental duplications on human chromosome 1. This gene family has experienced its greatest expansion within the human lineage and has expanded, to a lesser extent, among primates in general. Members of this gene family are characterized by tandemly repeated copies of DUF1220 protein domains. DUF1220 copy number variations in human chromosomal region 1q21.1, where most DUF1220 domains are located, have been implicated in a number of developmental and neurogenetic diseases such as microcephaly, macrocephaly, autism, schizophrenia, cognitive disability, congenital heart disease, neuroblastoma, and congenital kidney and urinary tract anomalies. Altered expression of some gene family members is associated with several types of cancer. This gene family contains numerous pseudogenes. [provided by RefSeq, Feb 2013]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <a href="mailto:techsupport@origene.com">techsupport@origene.com</a> . If you need a special design or shRNA sequence, please utilize our <a href="#">custom shRNA service</a> .



[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).