

## Product datasheet for **TL303023**

### Neuron navigator 1 (NAV1) Human shRNA Plasmid Kit (Locus ID 89796)

#### Product data:

Product Type:	shRNA Plasmids
Product Name:	Neuron navigator 1 (NAV1) Human shRNA Plasmid Kit (Locus ID 89796)
Locus ID:	89796
Synonyms:	POMFIL3; STEERIN1; UNC53H1
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	NAV1 - Human, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 89796). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	<a href="#">NM_001167738</a> , <a href="#">NM_020443</a> , <a href="#">NM_020443.1</a> , <a href="#">NM_020443.2</a> , <a href="#">NM_020443.3</a> , <a href="#">NM_020443.4</a> , <a href="#">NM_001167738.1</a> , <a href="#">BC007523</a> , <a href="#">BC069250</a> , <a href="#">BC117683</a> , <a href="#">BC156464</a> , <a href="#">BC172534</a> , <a href="#">NM_001167738.2</a>
UniProt ID:	<a href="#">Q8NEY1</a>
Summary:	This gene belongs to the neuron navigator family and is expressed predominantly in the nervous system. The encoded protein contains coiled-coil domains and a conserved AAA domain characteristic for ATPases associated with a variety of cellular activities. This gene is similar to unc-53, a Caenorhabditis elegans gene involved in axon guidance. The exact function of this gene is not known, but it is thought to play a role in neuronal development and regeneration. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Nov 2009]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <a href="mailto:techsupport@origene.com">techsupport@origene.com</a> . If you need a special design or shRNA sequence, please utilize our <a href="#">custom shRNA service</a> .



[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).