

Product datasheet for **TL303197**

KIF20B Human shRNA Plasmid Kit (Locus ID 9585)

Product data:

Product Type:	shRNA Plasmids
Product Name:	KIF20B Human shRNA Plasmid Kit (Locus ID 9585)
Locus ID:	9585
Synonyms:	CT90; KRMP1; MPHOSPH1; MPP-1; MPP1
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	KIF20B - Human, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 9585). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	NM_001284259 , NM_016195 , NM_016195.1 , NM_016195.2 , NM_016195.3 , NM_001284259.1 , BC005221 , BC012585 , BC042130 , BC046134 , BC058913 , BC093089 , BC108688 , NM_001284259.2 , NM_016195.4
UniProt ID:	Q96Q89
Summary:	Plus-end-directed motor enzyme that is required for completion of cytokinesis (PubMed:11470801, PubMed:12740395). Required for proper midbody organization and abscission in polarized cortical stem cells. Plays a role in the regulation of neuronal polarization by mediating the transport of specific cargos. Participates in the mobilization of SHTN1 and in the accumulation of PIP3 in the growth cone of primary hippocampal neurons in a tubulin and actin-dependent manner. In the developing telencephalon, cooperates with SHTN1 to promote both the transition from the multipolar to the bipolar stage and the radial migration of cortical neurons from the ventricular zone toward the superficial layer of the neocortex. Involved in cerebral cortex growth (By similarity). Acts as an oncogene for promoting bladder cancer cells proliferation, apoptosis inhibition and carcinogenic progression (PubMed:17409436).[UniProtKB/Swiss-Prot Function]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com . If you need a special design or shRNA sequence, please utilize our custom shRNA service .




[View online »](#)


**Performance
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.


For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).

Product images:


GFP signal was observed under microscope at 48 hours after transduction of TL303197A virus into HEK293 cells. TL303197A virus was prepared using lenti-shRNA TL303197A and [TR30037] packaging kit.



GFP signal was observed under microscope at 48 hours after transduction of TL303197B virus into HEK293 cells. TL303197B virus was prepared using lenti-shRNA TL303197B and [TR30037] packaging kit.



GFP signal was observed under microscope at 48 hours after transduction of [TL303197C] virus into HEK293 cells. [TL303197C] virus was prepared using lenti-shRNA [TL303197C] and [TR30037] packaging kit.



GFP signal was observed under microscope at 48 hours after transduction of [TL303197D] virus into HEK293 cells. [TL303197D] virus was prepared using lenti-shRNA [TL303197D] and [TR30037] packaging kit.