

Product datasheet for **TL303252V**

MITF Human shRNA Lentiviral Particle (Locus ID 4286)

Product data:

Product Type:	shRNA Lentiviral Particles
Product Name:	MITF Human shRNA Lentiviral Particle (Locus ID 4286)
Locus ID:	4286
Synonyms:	bHLHe32; CMM8; COMMAD; MI; WS2; WS2A
Vector:	pGFP-C-shLenti (TR30023)
Format:	Lentiviral particles
Components:	MITF - Human shRNA lentiviral particles (4 unique 29mer target-specific shRNA, 1 scramble control), 0.5 ml each, >10 ⁷ TU/ml.
RefSeq:	BC011461 , BC012503 , BC026961 , NM_000248 , NM_001184967 , NM_001184968 , NM_006722 , NM_198158 , NM_198159 , NM_198177 , NM_198178 , NM_001354604 , NM_001354605 , NM_001354606 , NM_001354607 , NM_001354608 , NM_198158.1 , NM_198158.2 , NM_000248.1 , NM_000248.2 , NM_000248.3 , NM_198159.1 , NM_198159.2 , NM_198177.1 , NM_198177.2 , NM_006722.1 , NM_006722.2 , NM_198178.1 , NM_198178.2 , NM_001184968.1 , NM_001184967.1 , BC026961.2 , BC065243 , BC065243.1 , BM800230 , NM_198178.3 , NM_001184967.2 , NM_000248.4 , NM_006722.3 , NM_001184968.2 , NM_198177.3 , NM_198158.3 , NM_198159.3
UniProt ID:	O75030
Summary:	The protein encoded by this gene is a transcription factor that contains both basic helix-loop-helix and leucine zipper structural features. The encoded protein regulates melanocyte development and is responsible for pigment cell-specific transcription of the melanogenesis enzyme genes. Heterozygous mutations in the this gene cause auditory-pigmentary syndromes, such as Waardenburg syndrome type 2 and Tietz syndrome. [provided by RefSeq, Aug 2017]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com . If you need a special design or shRNA sequence, please utilize our custom shRNA service .




[View online »](#)


**Performance
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).

Product images:

GFP signal was observed under microscope at 48 hours after transduction of TL303252A virus into HEK293 cells. TL303252A virus was prepared using lenti-shRNA TL303252A and [TR30037] packaging kit.



GFP signal was observed under microscope at 48 hours after transduction of TL303252B virus into HEK293 cells. TL303252B virus was prepared using lenti-shRNA TL303252B and [TR30037] packaging kit.