

Product datasheet for **TL303309**

MCM8 Human shRNA Plasmid Kit (Locus ID 84515)

Product data:

Product Type:	shRNA Plasmids
Product Name:	MCM8 Human shRNA Plasmid Kit (Locus ID 84515)
Locus ID:	84515
Synonyms:	C20orf154; dj967N21.5; POF10
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	MCM8 - Human, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 84515). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	NM_001281520 , NM_001281521 , NM_001281522 , NM_032485 , NM_182802 , NM_032485.1 , NM_032485.2 , NM_032485.3 , NM_032485.4 , NM_032485.5 , NM_182802.1 , NM_182802.2 , NM_001281522.1 , NM_001281520.1 , NM_001281521.1 , BC080656 , BC080656.1 , BC005170 , BC008830 , BC101054 , BC101055 , BC101056 , BC101057 , BM664368 , NM_182802.3 , NM_032485.6
UniProt ID:	Q9UJA3
Summary:	The protein encoded by this gene is one of the highly conserved mini-chromosome maintenance proteins (MCM) that are essential for the initiation of eukaryotic genome replication. The hexameric protein complex formed by the mini-chromosome maintenance proteins is a key component of the pre-replication complex and may be involved in the formation of replication forks and in the recruitment of other DNA replication related proteins. This protein contains the central domain that is conserved among the mini-chromosome maintenance proteins. The encoded protein may interact with other mini-chromosome maintenance proteins and play a role in DNA replication. This gene may be associated with length of reproductive lifespan and menopause. Alternatively spliced transcript variants encoding distinct isoforms have been described. [provided by RefSeq, Jul 2013]



[View online »](#)

- shRNA Design:** These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com. If you need a special design or shRNA sequence, please utilize our [custom shRNA service](#).
- Performance Guaranteed:** OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.
- For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).