

Product datasheet for **TL303325**

MBNL1 Human shRNA Plasmid Kit (Locus ID 4154)

Product data:

Product Type:	shRNA Plasmids
Product Name:	MBNL1 Human shRNA Plasmid Kit (Locus ID 4154)
Locus ID:	4154
Synonyms:	EXP; MBNL
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	MBNL1 - Human, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 4154). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	NM_001314057 , NM_021038 , NM_207292 , NM_207293 , NM_207294 , NM_207295 , NM_207296 , NM_207297 , NM_021038.1 , NM_021038.2 , NM_021038.3 , NM_021038.4 , NM_207295.1 , NM_207292.1 , NM_207292.2 , NM_207293.1 , NM_207294.1 , NM_207297.1 , NM_207296.1 , BC043493 , BC043493.1 , BC005296 , BC050535 , BM686432 , NM_001363870 , NM_207297.2 , NM_207292.3 , NM_021038.5 , NM_207293.2 , NM_207294.2 , NM_207296.2 , NM_207295.2
UniProt ID:	Q9NR56
Summary:	This gene encodes a member of the muscleblind protein family which was initially described in <i>Drosophila melanogaster</i> . The encoded protein is a C3H-type zinc finger protein that modulates alternative splicing of pre-mRNAs. Muscleblind proteins bind specifically to expanded dsCUG RNA but not to normal size CUG repeats and may thereby play a role in the pathophysiology of myotonic dystrophy. Mice lacking this gene exhibited muscle abnormalities and cataracts. Several alternatively spliced transcript variants have been described but the full-length natures of only some have been determined. The different isoforms are thought to have different binding specificities and/or splicing activities. [provided by RefSeq, Sep 2015]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com . If you need a special design or shRNA sequence, please utilize our custom shRNA service .



[View online »](#)

**Performance
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).