

Product datasheet for **TL303402**

LY6G5C Human shRNA Plasmid Kit (Locus ID 80741)

Product data:

Product Type:	shRNA Plasmids
Product Name:	LY6G5C Human shRNA Plasmid Kit (Locus ID 80741)
Locus ID:	80741
Synonyms:	C6orf20; G5c; G5c, NG33, C6orf20, LY6G5CA, LY6G5CB; LY6G5CA; LY6G5CB; lymphocyte antigen-6 G5C; lymphocyte antigen 6 complex, locus G5C; lymphocyte antigen 6 complex G5C; NG33
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	LY6G5C - Human, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 80741). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	NM_001002848 , NM_001002849 , NM_025262 , NM_025262.1 , NM_025262.2 , NM_025262.3 , NM_001002848.1 , NM_001002849.1 , BC140342 , BC141637
UniProt ID:	Q5SRR4
Summary:	LY6G5C belongs to a cluster of leukocyte antigen-6 (LY6) genes located in the major histocompatibility complex (MHC) class III region on chromosome 6. Members of the LY6 superfamily typically contain 70 to 80 amino acids, including 8 to 10 cysteines. Most LY6 proteins are attached to the cell surface by a glycosylphosphatidylinositol (GPI) anchor that is directly involved in signal transduction (Mallya et al., 2002 [PubMed 12079290]).[supplied by OMIM, Mar 2008]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com . If you need a special design or shRNA sequence, please utilize our custom shRNA service .



[View online »](#)

**Performance
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).