

## Product datasheet for **TL303490**

### LRDD (PIDD1) Human shRNA Plasmid Kit (Locus ID 55367)

#### Product data:

Product Type:	shRNA Plasmids
Product Name:	LRDD (PIDD1) Human shRNA Plasmid Kit (Locus ID 55367)
Locus ID:	55367
Synonyms:	DKFZp434D229; leucine-rich and death domain containing; leucine-rich repeats and death domain containing; leucine rich repeat and death domain containing protein; MGC16925; PIDD
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	PIDD1 - Human, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 55367). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	<a href="#">NM_018494</a> , <a href="#">NM_145886</a> , <a href="#">NM_145887</a> , <a href="#">NM_145887.1</a> , <a href="#">NM_145887.2</a> , <a href="#">NM_145887.3</a> , <a href="#">NM_145886.1</a> , <a href="#">NM_145886.2</a> , <a href="#">NM_145886.3</a> , <a href="#">NM_018494.3</a> , <a href="#">BC014904</a> , <a href="#">BC014904.1</a> , <a href="#">NM_145887.4</a> , <a href="#">NM_145886.4</a>
UniProt ID:	<a href="#">Q9HB75</a>
Summary:	The protein encoded by this gene contains a leucine-rich repeat and a death domain. This protein has been shown to interact with other death domain proteins, such as Fas (TNFRSF6)-associated via death domain (FADD) and MAP-kinase activating death domain-containing protein (MADD), and thus may function as an adaptor protein in cell death-related signaling processes. The expression of the mouse counterpart of this gene has been found to be positively regulated by the tumor suppressor p53 and to induce cell apoptosis in response to DNA damage, which suggests a role for this gene as an effector of p53-dependent apoptosis. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Aug 2010]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <a href="mailto:techsupport@origene.com">techsupport@origene.com</a> . If you need a special design or shRNA sequence, please utilize our <a href="#">custom shRNA service</a> .



[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).