

## Product datasheet for **TL303498**

### LOXHD1 Human shRNA Plasmid Kit (Locus ID 125336)

#### Product data:

Product Type:	shRNA Plasmids
Product Name:	LOXHD1 Human shRNA Plasmid Kit (Locus ID 125336)
Locus ID:	125336
Synonyms:	DFNB77; LH2D1
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	LOXHD1 - Human, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 125336). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	<a href="#">NM_001145472</a> , <a href="#">NM_001145473</a> , <a href="#">NM_001173129</a> , <a href="#">NM_001308013</a> , <a href="#">NM_144612</a> , <a href="#">NM_144612.1</a> , <a href="#">NM_144612.2</a> , <a href="#">NM_144612.3</a> , <a href="#">NM_144612.4</a> , <a href="#">NM_144612.5</a> , <a href="#">NM_144612.6</a> , <a href="#">NM_001145472.1</a> , <a href="#">NM_001145472.2</a> , <a href="#">NM_001145473.1</a> , <a href="#">NM_001145473.2</a> , <a href="#">BC041860</a> , <a href="#">BC047720</a> , <a href="#">NM_001145472.3</a>
UniProt ID:	<a href="#">Q8IVV2</a>
Summary:	This gene encodes a highly conserved protein consisting entirely of PLAT (polycystin/lipoxygenase/alpha-toxin) domains, thought to be involved in targeting proteins to the plasma membrane. Studies in mice show that this gene is expressed in the mechanosensory hair cells in the inner ear, and mutations in this gene lead to auditory defects, indicating that this gene is essential for normal hair cell function. Screening of human families segregating deafness identified a mutation in this gene which causes DFNB77, a progressive form of autosomal-recessive nonsyndromic hearing loss (ARNSHL). Alternatively spliced transcript variants encoding different isoforms have been noted for this gene. [provided by RefSeq, Mar 2010]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <a href="mailto:techsupport@origene.com">techsupport@origene.com</a> . If you need a special design or shRNA sequence, please utilize our <a href="#">custom shRNA service</a> .


[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).