

## Product datasheet for **TL303526**

### Lin28 (LIN28A) Human shRNA Plasmid Kit (Locus ID 79727)

#### Product data:

Product Type:	shRNA Plasmids
Product Name:	Lin28 (LIN28A) Human shRNA Plasmid Kit (Locus ID 79727)
Locus ID:	79727
Synonyms:	CSDD1; LIN-28; lin-28A; LIN28; ZCCHC1
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	LIN28A - Human, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 79727). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	<a href="#">NM_024674</a> , <a href="#">NM_024674.1</a> , <a href="#">NM_024674.2</a> , <a href="#">NM_024674.3</a> , <a href="#">NM_024674.4</a> , <a href="#">NM_024674.5</a> , <a href="#">BC028566</a> , <a href="#">BC028566.2</a> , <a href="#">NM_024674.6</a>
UniProt ID:	<a href="#">Q9H9Z2</a>
Summary:	This gene encodes a LIN-28 family RNA-binding protein that acts as a posttranscriptional regulator of genes involved in developmental timing and self-renewal in embryonic stem cells. The encoded protein functions through direct interaction with target mRNAs and by disrupting the maturation of certain miRNAs involved in embryonic development. This protein prevents the terminal processing of the LET7 family of microRNAs which are major regulators of cellular growth and differentiation. Aberrant expression of this gene is associated with cancer progression in multiple tissues. [provided by RefSeq, Sep 2015]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <a href="mailto:techsupport@origene.com">techsupport@origene.com</a> . If you need a special design or shRNA sequence, please utilize our <a href="#">custom shRNA service</a> .



[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).