

Product datasheet for **TL303625V**

Cytokeratin 6C (KRT6C) Human shRNA Lentiviral Particle (Locus ID 286887)

Product data:

Product Type:	shRNA Lentiviral Particles
Product Name:	Cytokeratin 6C (KRT6C) Human shRNA Lentiviral Particle (Locus ID 286887)
Locus ID:	286887
Synonyms:	K6E; KRT6E; PPKNEFD
Vector:	pGFP-C-shLenti (TR30023)
Format:	Lentiviral particles
Components:	KRT6C - Human shRNA lentiviral particles (4 unique 29mer target-specific shRNA, 1 scramble control), 0.5 ml each, >10 ⁷ TU/ml.
RefSeq:	NM_173086 , NM_173086.1 , NM_173086.2 , NM_173086.3 , NM_173086.4 , BC110639 , BC130583 , BC130585
UniProt ID:	P48668
Summary:	Keratins are intermediate filament proteins responsible for the structural integrity of epithelial cells and are subdivided into epithelial keratins and hair keratins. The type II keratins are clustered in a region of chromosome 12q13. [provided by RefSeq, Jul 2009]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com . If you need a special design or shRNA sequence, please utilize our custom shRNA service .



[View online »](#)

**Performance
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).