

Product datasheet for **TL303635V**

Cytokeratin 18 (KRT18) Human shRNA Lentiviral Particle (Locus ID 3875)

Product data:

Product Type:	shRNA Lentiviral Particles
Product Name:	Cytokeratin 18 (KRT18) Human shRNA Lentiviral Particle (Locus ID 3875)
Locus ID:	3875
Synonyms:	CK-18; CYK18; K18
Vector:	pGFP-C-shLenti (TR30023)
Format:	Lentiviral particles
Components:	KRT18 - Human shRNA lentiviral particles (4 unique 29mer target-specific shRNA, 1 scramble control), 0.5 ml each, $>10^7$ TU/ml.
RefSeq:	NM_000224 , NM_199187 , NM_199187.1 , NM_000224.1 , NM_000224.2 , BC072017 , BC072017.1 , BC000180 , BC000698 , BC004253 , BC008636 , BC009754 , BC020982 , NM_000224.3
UniProt ID:	P05783
Summary:	KRT18 encodes the type I intermediate filament chain keratin 18. Keratin 18, together with its filament partner keratin 8, are perhaps the most commonly found members of the intermediate filament gene family. They are expressed in single layer epithelial tissues of the body. Mutations in this gene have been linked to cryptogenic cirrhosis. Two transcript variants encoding the same protein have been found for this gene. [provided by RefSeq, Jul 2008]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com . If you need a special design or shRNA sequence, please utilize our custom shRNA service .


[View online »](#)

**Performance
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).